


Digital Body Weight Scale

User Manual

Read this manual thoroughly before using
and save it for future reference

Item Number: 58008

What's Included




1x Scale



2x AAA Battery

At a Glance



Setting Up Your Scale

- 1) Remove the plastic strip from the battery compartment.
- 2) Place the scale on a hard, flat surface.
- 3) Step on the scale and weigh yourself.

Note:

For best results, remember to place your scale on a hard, level surface. Step on the scale with bare feet and evenly distribute your weight.

Auto-on/off

To turn on the scale, simply step on the weighing platform with both feet and evenly distribute your weight. Once a steady reading is taken, the scale will lock the weight and the LCD will flash to indicate the final weight. .

The scale will automatically switch off in approximately eight (8) seconds without operation. After you have finished using the scale, please store the scale on a flat surface and not on its sides (refer to FAQ).

Unit Switch

Press the unit switch button on the bottom of the scale to alternate between Kilogram (kg), Pound (lb), and Stone (st).

Product Specs

Dimension:	11.8x11.8x0.7 inch / 300x300x18 mm
Weighing Platform:	8mm Tempered Glass
LCD Size:	3.1 inch / 74x28 mm
Display Type:	Backlit LCD
Battery:	2x AAA
Division:	0.1kg / 0.2lb
Weight Capacity: (Maximum)	180kg / 400lb / 28st

LCD Error Code

Lo

Low battery:
Replace the batteries before using again.

Err

Overload error:
Do not weight any object heavier than 400 lb.


E

Calculating error:
Step off and try again.

Warnings

- To avoid tipping, please be careful so as to not step on the edge of the scale.
- Do not jump on scale.
- Protect scale from hard hits, temperature fluctuations, and nearby heat sources (e.g. stoves, heating units).
- Do not drop the scale or other objects onto it as it may damage the internal sensors.
- Do not use the scale when your feet or body is wet.
- Do not use scale if it is wet.
- Do not place objects heavier than 400 lb/180 kg on the scale.
- Place the scale on level flooring-otherwise, it may tip.
- Clean the scale with a slightly damp cloth. DO NOT immerse the scale in water or use chemical/abrasive cleaning agents.

Battery Handling



- If the scale will be stored for a long period, it is recommended to remove the batteries to avoid leakage.

- Keep the batteries out of reach of children.

- Do not mix old and new batteries or batteries of different composition or brand.



- Do not heat or deform batteries.

- Do not dispose of batteries in normal waste.

- Check with your local authority for battery recycling advice.



FAQ

My scale sometimes shows a weight that is 4 to 5 pounds heavier than it should. It finally registers the correct weight after 2 to 3 tries.

Answer:

This can happen if you store your scale on its side, or for prolonged periods. Our scale has an "auto-on" feature which allows users to just step on and weigh without first calibrating the scale. The scale calibrates itself when it turns off, so if it is stored after use the scale's weight will not be properly accounted for, and may increase your reading by 4-5 pounds.

We understand that everyone doesn't have enough room in their bathroom to leave the scale out at all times. To account for this issue, the scale can be manually calibrated before weighing yourself by simply pressing the scale awake and waiting for it to zero out (0.0). Once zeroed out, the scale is ready for use.

Thank
You

Scan for more
product information



www.TenergyLife.com



Electrical waste should not be disposed with household waste.
Check with your local authority or retailer for recycling advice.